


Bottles for Injectables PP

100ml, 250ml, 500ml, 1000ml

PLASTIC
Compatible
with the
environment



| Ref. n° | 1894 | 1697B | 1698B | 1699B |
|---------------------------------|---------------------------------------------------------------------------------------------------|---------------|---------------|---------------|
| Volume (ml) | 100 | 250 | 500 | 1000 |
| Material | PP | PP | PP | PP |
| Colour | Natural/Amber | Natural/Amber | Natural/Amber | Natural/Amber |
| Dimension A (mm) | 98.5 | 121.8 | 157.6 | 203.3 |
| Dimension B (mm) | 64.4 | 77.5 | 106.2 | 141.8 |
| Dimension C (mm) | Ø49 | Ø66 | Ø78 | Ø95 |
| Neck opening (mm) | Ø22.5 | Ø22.5 | Ø22.5 | Ø22.5 |
| Weight (gr) | 17 | 30 | 43 | 62 |
| Brimful Volume | 127.0 | 279.0 | 531.0 | 1053.0 |
| Packing (pcs/brick) | 36 | 35 | 27 | 16 |
| Packing (pcs/carton) | 432 | 210 | 135 | 64 |
| Volume of carton (m3) | 0.144 | 0.144 | 0.144 | 0.144 |
| Carton dimension l x w x h (mm) | 800x400x450 | 800x400x450 | 800x400x450 | 800x400x450 |
| Carton weight filled (kg) | 8.3 | 7.4 | 7.3 | 5.5 |
| Cartons/EUR pallet | 15 | 15 | 15 | 15 |
| Quantity/EUR pallet | 6,480 | 3,150 | 2,025 | 960 |
| Sterilization | Autoclavable | | | |
| Rubber disc or stopper - Cap | The bottles fit rubber plugs according to ISO 8536-2 and aluminium caps according to ISO 8536-3/7 | | | |

2017-09-2

The data set forth herein has been carefully compiled by MedicoPack A/S. However, there is no warranty of any kind either express or implied, applicable to its use and the user assumes all risk and liability in connection therewith.