


Bottles for Injectables Polyethylene (HDPE)

250 mL, 500 mL



Ref. no.	1697	1698
Volume (mL)	250	500
Material	HDPE	HDPE
Colour	Natural	Natural
Dimension A (mm)	121.8	157.6
Dimension B (mm)	77.5	106.2
Dimension C (mm)	Ø66	Ø78
Neck opening (mm)	Ø22.5	Ø22.5
Weight (gr)	33.5	38
Brimful Volume (mL)	269	535
Packing (pcs/pack unit)	35	27
Packing (pcs/carton)	210	135
Volume of carton (m ³)	0.144	0.144
Carton dimension LxWxH (mm)	800x400x450	800x400x450
Carton weight filled (kg)	7.7	5.8
Cartons/EUR pallet	12	12
Quantity/EUR pallet	2,520	1,620
Sterilization	Irradiation	
Rubber disc or stopper - Cap	The bottles fit rubber plugs according to ISO 8536-2 and aluminium caps according to ISO 8536-3/7	

2024-09

The data set forth herein has been carefully compiled by MedicoPack A/S. However, there is no warranty of any kind either express or implied, applicable to its use and the user assumes all risk and liability in connection therewith.