

Sterile Plastic Vials

Substitution for glass



MedicoPack Sterile Plastic Vials

- Manufactured under clean room conditions ISO class 7.
- Documented sterile and supplied Ready-To-Use.
- Available in Polypropylene (PP) and COC (ClearVial™).
- Minimizes the risk of breakage and delamination.
- Drug Master Files available.

MedicoPack is a leading global supplier of Plastic Vials for the pharmaceutical industry. We are experts in semi-transparent Polypropylene (PP) vials and 100%-transparent COC vials (ClearVial™) – An alternative to glass vials to minimize the risk of breakage and delamination.



PP vials packed in sealed bags

Ref. no.	Volume (mL)	Diameter \varnothing (mm)	Height (mm)	Crimp finish \varnothing (mm)	Pcs per pack unit	Pcs per carton
1835ST01	5	22	37.9	20	230	4,140
1836ST01	10	26	45.6	20	162	2,592
1837ST01	20	32	48.6	20	98	1,372
1838ST01	25	33	53.5	20	85	1,105
1917ST01	50	43	71.5	20	55	940
1840ST01	100	49	87.6	20	41	492
1997ST01	250	66	121.8	32	35	210



PP vials packed in trays

Ref. no.	Volume (mL)	Diameter \varnothing (mm)	Height (mm)	Crimp finish \varnothing (mm)	Pcs per tray	Pcs per carton
1835ST51	5	26	37.9	20	228	2,280
1836ST51	10	26	45.6	20	160	1,600
1837ST51	20	32	48.6	20	104	936
1838ST51	25	33	53.3	20	104	832
1917ST51	50	43	71.5	20	54	378
1840ST51	100	49	87.6	20	40	400
1894ST51	100	49	98.5	32	40	360



COC vials packed in trays

Ref. no.	Volume (mL)	Diameter \varnothing (mm)	Height (mm)	Crimp finish \varnothing (mm)	Pcs per tray	Pcs per carton
1875T911	5	22	40	20	160	2,560
1876T911	10	24	45	20	149	2,384
1877T911	20	30	55	20	98	1,372
1879T911	50	42.5	73	20	50	500
1880T911	100	51.6	94.5	20	30	240